# **Corporate Information**

# **BOARD OF DIRECTORS**

Mr. Yao Yuan Mr. Iu Chung Mr. Chien Hoe Yong, Henry Mr. Hu Jun Mr. Yu Ti Jun Dr. Lam Lee G.\* Mr. Hu Jin Hua \* Mr. Lee Sze Ho, Henry \* \*Independent Non-Executive Directors

# AUDIT COMMITTEE

Dr. Lam Lee G. Mr. Hu Jin Hua Mr. Lee Sze Ho, Henry

#### AUTHORISED REPRESENTATIVES

Mr. Chien Hoe Yong, Henry Mr. Poon Kwong Wai, Kenny

#### **COMPANY SECRETARY**

Mr. Poon Kwong Wai, Kenny

## **AUDITORS**

Deloitte Touche Tohmatsu

# **LEGAL ADVISERS**

Freshfields Bruckhaus Deringer Baker & McKenzie Stephenson Harwood & Lo Chan & Tsu, Solicitors

#### **PRINCIPAL BANKERS**

DBS Bank (Hong Kong) Limited Citic Ka Wah Bank Limited The Hongkong and Shanghai Banking Corporation Limited Bank of China (Hong Kong) Limited Taishin International Bank Hong Kong Branch

#### PRINCIPAL SHARE REGISTRARS

Butterfield Corporate Services Limited 65 Front Street, Hamilton, Bermuda

#### **BRANCH SHARE REGISTRARS IN HONG KONG**

Central Registration Hong Kong Limited Shops 1712-1716, 17th Floor, Hopewell Centre, 183 Queen's Road East, Wanchai, Hong Kong

#### **REGISTERED OFFICE**

Canon's Court, 22 Victoria Street, Hamilton HM 12, Bermuda

#### **PRINCIPAL PLACE OF BUSINESS**

Room 1801-03, Hutchison House, 10 Harcourt Road, Central, Hong Kong Tel: (852) 3102 3201 Fax: (852) 3102 0905 Email: mingyuan@mingyuan-hk.com Website: www.mymedicare.com.hk

#### PLACE OF SHARE LISTING

The Stock Exchange of Hong Kong Limited

# **STOCK CODES**

The Stock Exchange of Hong Kong Limited: 233 Reuters: 233.HK Bloomberg: 233 HK

The Quality Control process of the C-12 protein chip: A technician is adding chemical substrates to the C-12 protein chip for detecting tumor markers effectiveness.

The HD-2001A ChipReader system developed by Shanghai HealthDigit. On the left is an optical system for the capture of signals from the protein chip. On the right is a computer system for data analysis.



# Mingyuan Medicare Development Company Limited

A microarray inkjet printer which is used for the fabrication of the C-12 protein chip. The system is being cleansed before arraying proteins. A lamp and a camera (upper left) are used to monitor the quality of the printing process.

